


Unsurpassed Ultra-Widefield image quality





## iCare EIDON Family

# Unsurpassed Ultra-Widefield image quality

iCare EIDON Family is a family of fundus imaging devices designed to provide Ultra-Widefield high resolution images. iCare EIDON Family is composed of three models, iCare EIDON, EIDON AF and EIDON FA, to serve different needs with a complete suite of imaging modalities.



TrueColor Confocal Technology, with white LED, promotes detailed retinal pictures with pupils as small as 2.5 mm and scanning through cataract.



**Ultra-Widefield\* high resolution imaging** provides both detail-rich and panoramic views of the retina up to 200°.



Flexibility in operational modes, from fully automatic control to manual operation, requires minimal staff training, speeding up the examination time.


# The TrueColor Confocal Technology

iCare EIDON Family combines a confocal optical engine with a white light LED source - which includes the entire visible spectrum of light - to illuminate the retina and to capture TrueColor fundus images characterized by colors close to reality.



# iCare EIDON Family allows a panoramic view of the retina, up to 200°


By using a proprietary stitching algorithm, iCare EIDON Family automatically combines different fields without user intervention, offering unsurpassed image quality, sharpness, and details, from the posterior pole to the periphery.




\*Ultra-Widefield imaging is available with the optional EIDON UWF Module.

## A complete suite of retina images


**TrueColor** to see the retina with superior color fidelity




Infrared to get information corresponding to the choroid



**Blue Fundus Autofluorescence** to assess the Retinal Pigment Epithelial layer integrity



**Fluorescein Angiography** to observe and monitor retinal blood flow details




All these imaging modalities are available also using the optional EIDON UWF Module, for wider field of view acquisition.


#### **RGB channels** separation

iCare EIDON Family TrueColor images can be separated into red, green and blue channel images to enhance the visualization of details on different retina layers.


TrueColor




Red Channel




Green Channel



Blue Channel



All these imaging modalities are available also using the optional EIDON UWF Module, for wider field of view acquisition.



### Stereo Viewer Technology


iCare EIDON Family Stereo Viewer Technology automatically captures two separate photos of the nasal field, at different angles and different focal planes (bifocal), creating improved 3D perception of the disc.

# **Fluorescein Angiography video** provides detailed dynamic view of retina vasculature

The fully automated capability of iCare EIDON FA transforms the overall complex Fluorescein Angiography process in a simple examination requiring minimal operator involvement.

Additionally, the Fluorescein Angiography video acquisition offers a dynamic perspective of retinal vasculature and related circulation mechanisms that may be missed with static photography.

With these functionalities, iCare EIDON FA leaves the operator free to focus attention on ensuring maximum patient care and comfort.



#### The new Ultra-Widefield Module


The early signs of several pathologies are often subtle and may appear first on the periphery of the retina.

The iCare EIDON UWF Module expands the retina field of view, revealig new areas that would not have been discovered with a standard field.




#### Preserved ultra-high image quality

Thanks to iCare EIDON UWF Module the iCare TrueColor Confocal Technology can be applied to Ultra-Widefield imaging, expanding the iCare EIDON Family's unsurpassed image quality to increased field of view on the retina.



### Increased field of view up to 200°

With iCare EIDON UWF Module it is possible to get a wider view of the fundus, imaging 120° of the retina in a single shot and even up to 200° with the Mosaic functionality, thus revealing signs of pathology outside the standard field and allowing for more thorough detection and documentation of peripheral retinal pathologies.



The 200° Infrared image shows clearly the vortex vein ampullae in all four retina quadrants.

### Provide your patient with a comfortable exam experience

iCare EIDON Family provides the flexibility of going from fully automatic to fully manual operation, requiring minimal staff training and speeding up examination time. With no need for dilation (2.5 mm minimum pupil), iCare EIDON Family promises enhanced comfort to patients, improving everyday workflow.



#### iCare EIDON Family Remote Viewer\*

The EIDON Family Remote Viewer is a browser-based software that allows for reviewing from any network computer on the same local area network (LAN), with password protection.



#### Technical data

	iCare EIDON	iCare EIDON AF	iCare EIDON FA
Fundus Imaging System Features			
Non-mydriatic Field of View (measured from the center of the eye) Individual image Mosaic*2 (three images) Sensor resolution	minimum pupil size 2.5 mm  standard 90° / Ultra-Widefield*¹ 120° standard up to 160° / Ultra-Widefield*¹ up to 200° 14 Mpixel (4608 x 3288)		
Light sources			
White LED (440 - 650 nm) Infrared LED (825 - 870 nm) Blue LED (440-475 nm)	<b>*</b>	* *	<b>*</b>
Imaging modalities			
TrueColor Infrared RGB channels separation* <sup>3</sup> Fundus Autofluorescence Fluorescein Angiography	<b>* * *</b>	* * * * * * * * * * * * * * * * * * *	* * * * * * * * * * * * * * * * * * *
Working distance Automatic operation Auto-focusing adjustment range Resolution Optical resolution on the retina Pixel pitch Fixation targets Dynamic Programmable Internal Fixation Target  FA video resolution	standard 28 mm / Ultra-Widefield* 16 mm auto-alignment, auto-focus, auto-exposure, auto-capture, auto-mosaic -12D to +15D 60 pixel/deg 15 microns 4.9 microns Internal/External  Central, Nasal, Temporal, Central-Nasal, Superior, Inferior, Superior-Temporal, Superior-Nasal, Inferior-Temporal, Inferior-Nasal		
FA video acquisition rate	5 fps		
Internal computer			
Hard drive Display Interfaces Export/Import		SSD, 480 GB nulti-touch, color display ort x 2, Gigabit Ethernet Port x	SSD, 2 TB
Dimensions/Power supply			
Size (WxHxD) Weight Related Voltage Frequency Power Consuption Electrical Class	360 mm X 590 mm X 6 25 kg (55 lbs) 100-240 VAC 50-60 Hz 80W IEC 60601-1 Class I	320 mm / 14.2" x 23.2" x 24.4"	

2021-

#### iCare. For better perception.

iCare is a trusted partner in ophthalmic diagnostics, offering physicians fast, easy-to-use, and reliable tools for diagnosis of glaucoma, diabetic retinopathy, and macular degeneration (AMD). Our product assortment includes automated TrueColor imaging devices, perimeters and handheld rebound tonometers.

We believe that ophthalmic care should be accessible, effortless, and reliable, and we aim to establish the next level of eye care.

EIDON, EIDON AF and EIDON FA are Class IIa Medical Devices, according to 93/42/EEC



#### Centervue S.p.A.

Via San Marco 9H 35129 Padova, Italy Ph. +39 049 501 8399 info@icare-world.com

#### **Icare Finland Oy**

Äyritie 22 01510 Vantaa, Finland Ph. +358 9 8775 1150 info@icare-world.com

#### Icare USA, Inc.

4700 Falls of Neuse Rd. Ste 245 Raleigh, NC. 27609 Ph. +1 888.422.7313 Fax +1 877.477.5485 infoUSA@icare-world.com

www.icare-world.com

